


Purchase Order No	Date	Fitter's Phone No.
Ordered by	Patient	Age
Street	Ship to	
City State Zip	Street	
Fitter	City State Zip	



Please refer to our current price list for a description of these items

**HOTLINE SERVICE**

(Guaranteed delivery in 4-5 business days, or less) **Additional Cost**

**REQUIRED Allergies:**

Allergic to silicone?  Yes,  No,  NA | Allergic to metal?  Yes,  No,  NA

**OPTIONAL SHOULDER STRAP MEASUREMENTS:** (See diagram)

**A. Length** from end of last strap to base of neck, with arm at side (same side as sleeve or gauntlet)

**B. Waist** Circumference

**C. Diagonal Length** for Shoulder Flap with retaining strap

**OPTIONS & EXTRAS:**

● **Zipper:**

- Lateral (Outside) (or)  Medial (Inside) (or)
- Posterior (Top)

● **Garment Color No:** \_\_\_\_\_

● **Elbow Lining:**  Inside (or)  Full (or)  Outside

● **Strap:**

- Worn to opposite axilla (or)
- Adjustable, axilla to waist (or)
- Two piece, with separate belt
- Bra strap, with Velcro flap

● **Flex/Contracture Seam:**  Elbow

● **Elastic Band:** Regular  1"  2" (or) Silicone  1"  2"

**DIAGNOSIS:** (Must be indicated to obtain proper counter-pressure)

- 22-28 mm Hg - Burns
- 25-30 mm Hg - Lymphedema Mild
- 30-35 mm Hg - Lymphedema Moderate
- 35-40 mm Hg - Lymphedema Severe

**COMMENTS & INSTRUCTIONS:**

---



---

**Check this box if there are instructions on other side**

Left	Right
3	
1½	
0	
1½	
3	
4½	
6	
7½	
9	
10½	
12	
13½	
15	
16½	
18	
19½	

**For Gauntlets Only**

\_\_\_\_\_ Circumference between thumb joints  
 Separate (or)  One Piece with Sleeve

**Wrist** - Measurement is taken ahead of wrist bone, toward hand

**Optional Shoulder Strap Measurements (A, B and C)**

**Elbow** - Anchor Point

