


Purchase Order No _____	Date _____	Fitter's Phone No. _____
Ordered by _____	Patient _____	Age _____
Street _____	Ship to _____	
City _____ State _____ Zip _____	Street _____	
Fitter _____	City _____ State _____ Zip _____	




PLEASE CHECK ONE: Burn | Lymphedema

For Body Suits #418 - #422: Complete girth and length measurements desired on Form F-029 or Form F-030.

For Arm Sleeves on Body Suits #412 - #422: Complete Form F-035 or Form F-036.

An oversize charge will be made for any girth measurement exceeding 50 inches.

Form F-040 - Body Measuring



Please refer to our current price list for a description of these items

HOTLINE SERVICE

(Guaranteed delivery in 4-5 business days, or less) **Additional Cost**

REQUIRED Allergies:

Allergic to silicone? Yes, No, NA | Allergic to metal? Yes, No, NA

Numbers - Measure Full Circumferences (All Dotted Lines)

- 1. Neck
- 2. Right shoulder sleeve seam
- 3. Left shoulder sleeve seam
- 4. Chest at axilla
- 5. Female support only - Chest at nipple line and cup size _____
- 6. Waist (standard)
- 7. Below waist (optional)
- 11. 2" below axilla (short sleeve only)

Following Girths for Catalog Numbers: 414 - 422

- 8. Largest part of buttocks
- 9. Leg at fold of buttocks

Letters - Length (Solid lines)

- A. Neck to shoulder tip
- B. (Female support only) base of neckline to under breasts (at sternum)
- C. Neck line to waist in front & back
- D. Neck line to below waist (if needed) in front & back

Following Lengths for Catalog Numbers: 414 - 417 (briefs only)

- E. Neck line to largest part of buttocks
- F. Neck to fold of buttocks Crotch: Snaps (or) Velcro

Following Lengths for Catalog Numbers: 418 - 422

- E. Neck line to largest part of buttocks
- F. Neck to fold of buttocks Crotch: Open (or) Closed

Options and Extras

- Arm Zipper, indicate length & location on form F-035 or F-036
- Leg Zipper, indicate length & location on form F-029 or F-030

● **Garment Color No:** _____

● **Axilla (Closed is standard)**

- Open Left Open Right Mesh Left Mesh Right
- Crewneck V - Neck Scoopneck Turtle Neck
- Lining inside elbow Lining full elbow Flexseam elbow
- Front Closure Zipper Velcro
- Back Closure Zipper Velcro

Check this box if there are instructions on other side