


Purchase Order No	Date	Fitter's Phone No.
Ordered by	Patient	Age
Street	Ship to	
City State Zip	Street	
Fitter	City State Zip	

**Letters (Sleeves & Gauntlets)**

- A. 3" from wrist (O)
- B. 6" from wrist (O)
- C. Elbow
- D. Halfway between C and E
- E. Axilla
- F. Palm at web of thumb
- G. Palm at base of fingers
- H. Between thumb joints
- O. Wrist, ahead of wrist bone toward hand




**SLEEVES**


LEFT			RIGHT		
Girth	Length		Girth	Length	
E	O to E		E	O to E	
D	O to D		D	O to D	
C	O to C		C	O to C	
B	O to B	6" (15.2 cm)	B	O to B	6" (15.2 cm)
A	O to A	3" (7.6 cm)	A	O to A	3" (7.6 cm)
O	X		O	X	

**GAUNTLETS**  Separate  One Piece with Sleeve

LEFT			RIGHT		
Girth	Length		Girth	Length	
F	O to F		F	O to F	
G	O to G		G	O to G	
H	X		H	X	
O	X		O	X	



**Optional Shoulder Strap Measurements**  
(1, 2 and 3)



Please refer to our current price list for a description of these items

**HOTLINE SERVICE**

(Guaranteed delivery in 4-5 business days, or less) **Additional Cost**

**REQUIRED Allergies:**

Allergic to silicone?  Yes,  No,  NA | Allergic to metal?  Yes,  No,  NA

**OPTIONAL SHOULDER STRAP MEASUREMENTS:** (See diagram)

1. **Length** from end of last measurement to base of neck, with arm at side (same side as sleeve or gauntlet)
2. **Waist Circumference**
3. **Diagonal Length** for Shoulder Flap with retaining strap

**OPTIONS & EXTRAS:**

- **Zipper:**  Lateral (Outside) (or)  Medial (Inside) (or)  Posterior (Top)
- **Garment Color No:** \_\_\_\_\_
- **Elbow Lining:**  Inside (or)  Full (or)  Outside
- **Strap:**
  - Worn to opposite axilla (or)
  - Adjustable, axilla to waist (or)
  - Two piece, with separate belt
- **Flex/Contracture Seam:**  Elbow
- **Elastic Band:** Regular  1"  2" (or) Silicone  1"  2" (or) Microdot  2"

**DIAGNOSIS:** (Must be indicated to obtain proper counter-pressure)

- 22-28 mm Hg - Burns
- 25-30 mm Hg - Lymphedema Mild
- 30-35 mm Hg - Lymphedema Moderate
- 35-40 mm Hg - Lymphedema Severe

**Check this box if there are instructions on other side**