


Purchase Order No	Date	Fitter's Phone No.
Ordered by	Patient	Age
Street	Ship to	
City State Zip	Street	
Fitter	City State Zip	

Letters (Sleeves & Gauntlets)

- A. 3" from wrist (O)
- B. 6" from wrist (O)
- C. Elbow
- D. Halfway between C and E
- E. Axilla
- F. Palm at web of thumb
- G. Palm at base of fingers
- H. Between thumb joints
- O. Wrist, ahead of wrist bone toward hand




SLEEVES

LEFT			RIGHT		
Girth	Length		Girth	Length	
E	O to E		E	O to E	
D	O to D		D	O to D	
C	O to C		C	O to C	
B	O to B	6" (15.2 cm)	B	O to B	6" (15.2 cm)
A	O to A	3" (7.6 cm)	A	O to A	3" (7.6 cm)
O	X		O	X	


GAUNTLETS Separate One Piece with Sleeve

LEFT			RIGHT		
Girth	Length		Girth	Length	
F	O to F		F	O to F	
G	O to G		G	O to G	
H	X		H	X	
O	X		O	X	



Optional Shoulder Strap Measurements
(1, 2 and 3)

Form F-036 IMF - for Use Without Paper Tapes



Please refer to our current price list for a description of these items

HOTLINE SERVICE

(Guaranteed delivery in 4-5 business days, or less) **Additional Cost**

REQUIRED Allergies:

Allergic to silicone? Yes, No, NA | Allergic to metal? Yes, No, NA

OPTIONAL SHOULDER STRAP MEASUREMENTS: (See diagram)

1. **Length** from end of last measurement to base of neck, with arm at side (same side as sleeve or gauntlet)
2. **Waist Circumference**
3. **Diagonal Length** for Shoulder Flap with retaining strap

OPTIONS & EXTRAS:

- **Zipper:** Lateral (Outside) (or) Medial (Inside) (or) Posterior (Top)
- **Garment Color No:** _____
- **Elbow Lining:** Inside (or) Full (or) Outside
- **Strap:**
 - Worn to opposite axilla (or)
 - Adjustable, axilla to waist (or)
 - Two piece, with separate belt
- **Flex/Contracture Seam:** Elbow
- **Elastic Band:** Regular 1" 2" (or) Silicone 1" 2"

DIAGNOSIS: (Must be indicated to obtain proper counter-pressure)

- 22-28 mm Hg - Burns
- 25-30 mm Hg - Lymphedema Mild
- 30-35 mm Hg - Lymphedema Moderate
- 35-40 mm Hg - Lymphedema Severe

Check this box if there are instructions on other side